



# Spengler

DEPUIS 1907



# Mobi

**PROFESSIONALS  
SPHYGMOMANOMETER  
& CUFF**



## Spengler

Spengler S.A.S | ZI La Limoise | Rue Robert Maréchal | 36100 ISSOUDUN | FRANCE

[www.spengler.fr](http://www.spengler.fr)

Tél. : +33 (0)4 42 52 89 80

Mail : [export@spengler.fr](mailto:export@spengler.fr)

# PROFESSIONAL SPHYGMOMANOMETER & CUFF

# Mobi



"LIFE LINK®" TUBE

The new "LIFE LINK®" tubing allows to change cuff in less than a second thanks to its magnetic connectors. Choose the cuff fitted to the patient's morphology and increase the accuracy



AMBIDEXTROUS

Specially designed to be ambidextrous, choose your configuration and change the position of the decompression valve, the connector "LIFE LINK®" and the spoon



ERGO PUMP

The ergonomic bulb is suitable for all types of hands and offers an extra volume of 11% compared to traditional bulbs



PRECISION VALVE

MOBI is equipped with a decompression valve made of a high sensibility metal alloy designed to provide high precision deflation



SHOCK RESISTANT

Due to its ABS plastic construction, MOBI is specially designed to offer great resistance to shocks and falls from everyday life



EASY CUFF

Developed in TPU coating, the MOBI cuff offers weight gain and easiness of cleaning



## DENOMINATION

- **Mobi Black Carbon with Adult Cuff (M)**
- **Mobi Blueberry with Adult Cuff (M)**
- **Mobi Coral red with Adult Cuff (M)**
- **Mobi Powdery pink with Adult Cuff (M)**

## REFERENCE

513 702  
513 700  
513 701  
513 703

## SPENGLER QUALITY

### CE MARKING

European directive 93/42/CEE  
mechanical sphygmomanometer : class 1 measuring function

Spengler sphygmomanometers meet the essential requirements of safety health, hygiene and environmental protection

### After sales service

Spengler provides maintenance and repair of equipment

### Legal Guarantee

Spengler sphygmomanometers are guaranteed for life against defects in workmanship hidden from the date of purchase

### Spengler guarantee

Spengler warranty of 3 years (under normal conditions of use): repair or replacement free parts recognized defective (malfunction due to a defect in material, workmanship or design) by our after sales service

### METROLOGICAL CONTROL

This blood pressure meter is a professional precision medical device. It is recommended to check the correct operation and accuracy every 2 years



[www.spengler.fr](http://www.spengler.fr)

Spengler S.A.S | ZI La Limoise | Rue Robert Maréchal | 36100 ISSOUDUN | FRANCE

Tél. : +33 (0)4 42 52 89 80

Mail : [export@spengler.fr](mailto:export@spengler.fr)